# SURGICAL INSTRUMENTS

**INSTRUCCIONES DE USO**  
**INSTRUCTIONS FOR USE**

**CONTENIDO**  
**NO ESTÉRIL**

**CONTENTS**  
**NOT STERILE**

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1. ANSI/AAMI ST79:2006 "Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities".

2. AAMI TIR12:2010 "Designing, testing, and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers."

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**IMPORTANT**

Please read the Instructions For Use carefully. The surgical instruments should not be utilized if the Instructions For Use have not been thoroughly read. Not following the indications of the Instructions For Use carefully may lead to the malfunction of the instruments and/or cause harm to the patient.

**INDICATIONS OF USE**

This surgical instrument manufactured by DIMA is intended to be used along with implantable products designed and manufactured by DIMA.

*Note: This product is provided non-sterile. New and used instruments must be thoroughly processed according to these instructions prior to use in order to attain sterility.*

**WARNINGS AND PRECAUTIONS**

- Surgical instruments should only be used by medical specialists trained in the use of these mechanisms. Surgeons should not begin clinical use of the surgical instrument without reviewing the Instructions for Use. It is imperative that appropriate personnel read all guideline materials prior to surgery. If the guidelines are not properly followed, the risk of complications to the patient increases.
- Before using the product, the packaging must be inspected to guarantee that it has not suffered any damages.
- Universal precautions should be observed by all hospital personnel that might be working with the contaminated or potentially contaminated surgical devices. Personal protective equipment (PPE) should be worn when handling or working with contaminated or potentially contaminated medical devices and equipment. PPE should include, among other options: gown, facial mask, goggles or face shield, gloves and shoe covers.
- Metal brushes or scouring pads must not be used during manual cleaning procedures. These materials will damage the surface and finish of instruments. Soft-bristled, nylon brushes and pipe cleaners should be used.
- Cleaning agents with low foaming surfactants should be used during manual cleaning procedures to ensure that instruments are visible in the cleaning solution. Manual scrubbing with brushes should always be performed with the instrument below the surface of the cleaning solution to prevent formation of aerosols and splashing which may spread contaminants. Cleaning agents must be easily and completely rinsed from device surfaces to prevent accumulation of detergent residue.
- Do not allow contaminated devices to dry prior to reprocessing. All subsequent cleaning and sterilization steps are facilitated by not allowing blood, body fluid, tissue debris, saline, or disinfectants to dry on used instruments.
- The instrument must be manually cleaned immediately after use. This initial cleaning should not be postponed until the cleaning and resterilization procedure.
- Saline and cleaning/disinfecting agents containing aldehyde, mercury, active chlorine, bromine, bromide, iodine or iodide are corrosive and should not be used; instruments must not be placed or soaked in Ringers Solution
- Descaling agents that include morpholine should not be used in steam sterilizers. These agents leave residue which can damage polymer instruments over time.
- As with any surgical instrument, careful attention should be made to assure that excessive force is not placed on this instrument. Excessive force can result in failure.
- If the instrument does not function properly, immediately discontinue use.
- Follow instructions and warnings as issued by manufacturers of any decontaminants, disinfectants and cleaning agents used. Wherever possible, avoid the use of mineral acids and harsh, abrasive agents.
- When reprocessing medical devices, always handle with care and wear protective clothing, gloves and eyewear in accordance with applicable health and safety procedures.
- No part of the process shall exceed 145°C
- The date of the last revision of this document is included at the end of the Instructions For Use. If more than 24 months have passed since the last revision, please contact the manufacturer or distributor to receive the updated Instructions for Use.

**MATERIALS USED IN THE MANUFACTURE OF THE INSTRUMENT(S)**

- Stainless steel
- Polypropylene (in the case that the instrument(s) contain non-metal parts)

**LIMITATIONS AND RESTRICTIONS ON REPROCESSING**

Repeated processing has minimal effect on these instruments. End of life is normally determined by wear and damage due to use.

**CARE AND MAINTENANCE INSTRUCTIONS**

- Proper care and maintenance is important for the efficient and safe operation of the medical/surgical equipment. Careful inspection of all equipment is highly recommended upon receiving it and prior to each use, as a safeguard against possible injury to patient and/or medical personnel operating the device or equipment.
- Visually inspect the instrument for damage or wear prior to and after cleaning. Discard instruments with damage or excessive wear.
- The distal tip and slot features should be smooth and without denting, cracks, burrs or fractured surfaces.
- Any instrument including a handle should be thoroughly examined to assure that needle is securely connected to the handle. All surfaces and edges should be uniform and smooth, free of corrosion and dents, sharp edges or points. Minor surface hazing and scratches on the wire are acceptable.

**PREPARATION FOR USE**

Instrument(s) provided non-sterile.

Before first and following uses, the instrument(s) must be cleaned and sterilized. Sterilization is done by steam autoclave. To ensure the usability and reliability of the product, it must be cleaned after every surgery. Both a manual and an automatic cleaning method are proposed in these instructions for use. Instrument(s) can be reused after conducting appropriate cleaning and sterilization.

**CLEANING INSTRUCTIONS**

To prevent the drying of blood and tissue residues and growth of microorganisms, instruments should be immersed, as soon as possible after each procedure, in an aqueous solution and preferably a combined neutral or enzymatic cleaning and/or disinfectant solution prepared and maintained per instructions of the solution manufacturer. This immersion process should be carried out as soon as reasonably practical following use of instrument(s). Regardless, the cleaning and resterilization process of the instrument(s) should also be performed as soon as reasonably possible, allowing for the reutilization of the instrument(s).

The instrument(s) sterilization is to be done by steam autoclave. To ensure the usability and reliability of the instrument(s), cleaning of equipment must be carried out after every intervention. Hereunder both a manual and an automatic cleaning method are proposed.

**MANUAL CLEANING PROCEDURE**

Upon termination of use, thoroughly rinse instrument(s) under cool running tap water (~20°C) for at least 90 seconds to remove all organic residue and contaminants. Use of a non-metallic soft bristle brush for light scrubbing is acceptable, paying special attention to areas where contaminants may be trapped (corners, grooves and holes.)

1. For 10 minutes, fully immerse the instrument(s) in enzymatic detergent solution for stainless steel (following the instructions for use of the detergent manufacturer) using lukewarm (approximately ~33°C) tap water. Verify that the expiration date of the detergent has not passed.
2. Thoroughly clean the instrument, scrubbing with a non-metallic bristle brush. This part of the cleaning process should be executed with the instrument(s) fully

- submerged in the enzymatic solution to avoid splashing, paying special attention to areas where contaminants may be trapped (corners, grooves, crevices and/or lumens).
3. Rinse the instrument(s) thoroughly in cold (approximately 20°C), running tap water.
  4. For at least 10 minutes, fully immerse the instrument(s) in a disinfectant detergent solution (following the instructions for use of the detergent manufacturer) and lukewarm (approximately 33°C) running tap water.
  5. Thoroughly clean the instrument, scrubbing with a non-metallic bristle brush. This part of the cleaning process should be executed with the instrument(s) fully submerged in the detergent solution to avoid splashing, paying special attention to areas where contaminants may be trapped (corners, grooves, crevices and/or lumens).
  6. Rinse the instrument(s) thoroughly with an abundant amount of distilled water, paying special attention to the hard-to-reach areas where contaminants may remain (corners, grooves, crevices and/or lumens).
  7. Visually inspect the instrument for cleanliness, paying special attention to the hard-to-reach areas where contaminants may remain (corners, grooves, crevices and/or lumens).  
If any residue is detected, return to STEP 1 and repeat the cleaning process completely. Once repeated, if no residue is detected, continue on with the drying procedure.
  8. Instrument(s) must be dried as soon as possible after cleaning in a protected or clean environment to minimize new contamination. Use a cotton cloth to dry instrument(s).

#### **AUTOMATIC CLEANING PROCEDURE**

Instrument(s) may be cleaned using ultrasounds. Hereafter a possible cleaning cycle is described:

1. Insert the instrument(s) into the ultrasonic washing machine.
2. Fill the tank of the ultrasonic machine with detergent and running tap water with a temperature of 30°C to 35°C (86°F to 95°F). Initiate the ultrasounds and leave on for 10 minutes. Empty the tank of the ultrasonic machine.
3. Fill the tank of the ultrasonic machine with fresh detergent and running tap water with a temperature of 30°C to 35°C (86°F to 95°F). Initiate the ultrasounds and leave on for 10 minutes. Empty the tank of the ultrasonic machine.
4. Fill the tank of the ultrasonic machine with distilled water at a temperature between 30°C to 35°C (86°F to 95°F). Initiate the ultrasounds and leave on for 10 minutes. Empty the tank.
5. Fill the tank of the ultrasonic machine with fresh distilled water at a temperature between 30°C to 35°C (86°F to 95°F). Initiate the ultrasounds and leave on for 10 minutes. Empty the tank of the ultrasonic machine.
6. Dry at 93°C (199°F) for 30 minutes.

#### **STERILIZATION PROCEDURE**

Instrument(s) is proved non-sterile.

Before first and following uses, the instrument(s) must be cleaned and sterilized.

Instrument(s) can be reused after conducting appropriate cleaning and sterilization.

1. Instrument(s) will be properly packaged in order to guarantee the durability of the sterilization procedure.
2. Once instrument(s) are correctly packaged, they will be inserted into an autoclave to proceed with steam sterilization.
3. The person responsible for carrying out the sterilization process should carefully follow the instructions for use provided by the autoclave manufacturer, given that adequate times and temperatures to guarantee complete sterilization of instrument(s) can vary and depend on the indications outlined by each specific manufacturer.
4. Once the autoclave sterilization process has finished, the packaged instrument(s) will be removed from the autoclave and will need to cool down.
5. It is absolutely necessary to verify that the packaging of the sterilized instrument(s) is completely dry. If condensation is observed in or on the packaging, they are considered unsterile and must be reprocessed (sterilized).
6. Once the instrument(s) cool, they will then proceed to be stored in a place that guarantees the complete preservation of the packing and the instrument(s).

The operator should follow the instructions for steam sterilization provided by the sterilizer's manufacturer which indicate cycle parameters. The following sterilization times and temperatures have been validated for Dynamic air removal/Pre-vacuum Sterilization:

- 3 preconditioning pulses
- Exposure Temperature: 132°C
- Exposure Time: 5 minutes
- Drying Time: 30 minutes

The previously described procedure should be considered as a general recommendation given that the sterilization procedure should be meticulously established following the instructions provided by the autoclave manufacturer. The responsibility of developing a proper sterilization procedure for the instrument(s) lies with the organization which establishes the procedure, NOT with the manufacturer of the instrument(s).

#### **STORAGE**

It is recommended to store in a dry place at room temperature, away from heat and direct light.

#### **ADDITIONAL INFORMATION**

Precautions with respect to the transmission of Creutzfeldt-Jacob disease (CJD)

- CJD is the human transmissible spongiform encephalopathy (TSE) variant agent which is recognized as resistant to the normal disinfection and sterilization procedures.
- The intended use for instrument(s) is recognized as low risk procedures for the transmission of CJD. Deactivation of the CJD agent is being actively researched at this time and is beyond the scope of these instructions
- Facilities reprocessing surgical instruments must take particular precautions when handling instruments used on patients at high risk for CJD to not only deactivate the CJD agent but also to reduce cross-contamination during the transportation and cleaning processes.

DIMA SL as manufacturer and distributor of the instruments, which are the subject of this manual, is not responsible for any direct or indirect damages resulting from improper or incorrect use, care, or servicing of these instruments, or from the failure to follow the guidelines contained in this manual.

If non-authorized persons perform any repairs on or modifications of this product, DIMA SL does not accept any liability and the warranty becomes void. Components influencing safety may only be replaced with original replacement parts.

This product and manual are subject to technical changes and modifications without notice.

#### **References**

These recommendations are developed using the following guides:

1. ANSI/AAMI ST79:2006 "Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities".
2. AAMI TIR12:2010 "Designing, testing, and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers."  
*Note: All cleaning and sterilization procedures should be validated by the health care facility performing these operations.*