

# From JUNE: Fascia Needle Passers & Handles

## KNOWLEDGE AND SKILL

We have **been chosen** to supply over 50 trusts in the UK, and we hope we can be of service to you as well. **We will do our best.**

## ONLY WOMEN'S HEALTH

We are the **first and only** company in the UK whose sole focus is on women's health. It is all we do, and nothing distracts us.

## CLOSE AND CAPABLE

We are **appreciated** as a supplier because we have staff and emergency stock in 6 locations throughout the country. We have prevented the cancellation of many surgeries by hand delivering product in emergency situations.

## ENVIRONMENT

We are **mindful about our impact**, and will try to travel smart. We reuse our shipping boxes (our Ugly and Proud campaign!) and we do what we can to preserve our planet whilst fulfilling our obligations to our customers. Our full policy is published on our website.



PAS22 and handle MPAS11

### Quality in design

Strong, solid quality steel makes for a dependable needle passer and handle in these reusable instruments for fascia repairs. Designed to fit into the procedure without changing needle trajectory or curvature, the needles come in two different diameters: 3.5mm thin or a thicker 5mm. Surgeon preference determines the choice.

### Variety of patients

For patients with a higher BMI, it can sometimes be challenging to keep control and visibility of the needle as the layers are too thick: we have solved this by making the needles in two different lengths: 230mm and 280mm.

- **PATIENTS:** Designed for native tissue repair using autologous fascia to treat female incontinence.
- **ENVIRONMENT:** Reusable. Manufactured in Spain. CE marked. Single packed.
- **OPTIONS:** Both thin and thick needles available: 3.5mm or 5mm diameter needle passers. 2 different style of handle.
- **SAFETY:** For normal as well as higher BMI patients: 230mm or 280mm length.
- **READ MORE:** Images and videos on [www.JUNEMEDICAL.CO.UK/products](http://www.JUNEMEDICAL.CO.UK/products)

### Importance of the handle

There are strong voices in support for the T-style handle, arguing that it sits better in the hand and gives a better balanced feel for the trajectory of the needle tip.

### Your choice

However, those who prefer the slim line handle will point out that they feel it to be an elongation of the needle into their hand, giving them a better sense for where they are in the anatomy. Whichever one is your favourite, we have one that you will like.

Solid manufacturing, no product gapping or clicking, you will feel in complete control of the procedure.



<- PAS22 needle tip with opening for suture

## TECHNICAL SUPPORT

Our team is extremely **knowledgeable**: we employ engineers, midwives, nurses and scientists: we are always on hand to answer any questions.

## ZERO WAIT

Sometimes your equipment needs a **repair**: we will get you a replacement machine and fix yours in no-time. Need an answer to a **question**? We will get back to you within 24hrs.

## CHARITY

Our focus for charity is curing obstetric fistula. Read more here: [www.fistulafoundation.com](http://www.fistulafoundation.com) and on our website.

## MADE IN BRITAIN

Proudly British, we are headquartered in the South East, close to major airports and motorways. We strive to **innovate and produce** in the UK and we work closely with the NHS, KTN, ASHN and DIT to bring only the best medical devices to our customers.

For more information on any of our products or policies please visit [www.junemedical.co.uk](http://www.junemedical.co.uk)

## Further details



MPAS09 T-style handle

### How does it work?

Needles are connected to the handle with a simple and sturdy screw function to give a secure snug connection between handle and needle. Handles and needles are reusable and resterilisation instructions are in the pack. The needle has a hole for the suture to pass through and connect to the fascia sling.

### What is used?

Most commonly we see 2 needles and 1 handle on each tray, but we have some customers who prefer to have both long and short needles available. This ensures they always have the best length depending on their patients anatomy.

### How many do we need?

Usually three complete sets are purchased, allowing for up to three patients to be treated on a list.

\* image is showing T-style handle AR-05, 2 needle passers 3.5mm PAS22 connected to a slimline handle MPAS11

## Order codes & cost

Product Name	Details	Product Code	Price excl VAT
Passer <b>Needle</b> 3.5 mm 230 mm	Connects with handle MPAS09. Non sterile, reusable.	PAS21	£100
Passer <b>Needle</b> 3.5mm/280 mm	Connects with handle MPAS09. Non sterile, reusable.	PAS22	£110
<b>Handle</b> 3.5 mm (slimline)	For use with either PAS21 or PAS22 Non sterile, reusable.	MPAS11	£595
<b>Handle</b> 3.5 mm (T-style)	For use with either PAS21 or PAS22. Non sterile, reusable.	MPAS09	£490
Passer <b>Needle</b> 5mm/230 mm	Fits handle AR-05. Non sterile, reusable.	AR-03	£120
Passer <b>Needle</b> 5mm/280 mm	Fits handle AR-05. Non sterile, reusable.	PAS05	£130
<b>Handle</b> 5mm (T-style)	For use with AR-03 or PAS05 Non sterile, reusable.	AR-05	£595
<b>Suprapubic</b> 3.5mm (T-style)	For top-down use, handle and needle preconnected. Non sterile, reusable.	APAS401	£334

## ORDERING INFO

Order@junemedical.co.uk

## GET IN TOUCH

Good old phonecall:

01628-330010

Twitter: @junemedical

Email: [contact@junemedical.co.uk](mailto:contact@junemedical.co.uk)

## PRODUCT SUPPORT

(8am-5.30pm)

WhatsApp or FaceTime:

07898-780418



INNOV8, Queen Alexandra Road  
High Wycombe, Buckinghamshire  
HP11 2GZ, United Kingdom